



PRODUCT FEATURES:

- fast and accurate measurements of focal spot size

RaySafe Pro-Slit

The RaySafe Pro-Slit phantom is a so called slit camera for accurate measurement of the focal spot size according to IEC 60336:2005. Its smart design enables repeatability and accuracy and the possibility to measure the size of any focal spot with one tool.

The phantom fits perfect with the Pro-Stand (not included) in order to enable easy and repeatable measurement setup.

TECHNICAL DATA

- Dimensions: 40 x 40 x 11.5 mm
- Material: Made out of tungsten which is embedded in PMMA
- Slit size: 0.01 mm wide and 10 mm long slit with 8° spread

